

www.awaretech.com

LumiStat[®] 4900

Programmable Open System

CHEMILUMINESCENT Microstrip Reader



sales@awaretech.com

LumiStat[®], part of the Stat Fax[®] line of semi-automated instruments from Awareness Technology, Inc. is a CLIA microstrip reader capable of processing CLIA assays

The sensitive Luminometer with a broad dynamic range is designed to read luminescent reactions taking place in up to 3 x 8 or 3 x 12 CLIA microwell strips. LumiStat[®] is also capable of reading 1 or 2 wells if desired.

CLIA Assays can be programmed using a variety of calculation modes as well as an RLU/absorbance mode. All instrument adjustments and settings are made using the built-in touchscreen interface.

The photomultiplier light detector is sensitive to light in the UV and visible spectrum from 340 nm to 650 nm.

LumiStat[®] is a stand-alone analyzer with the ability to transmit and capture data to a PC via USB connectivity.

Features include an incubation timer, automatic blanking options, the ability to indicate the locations of positive and negative controls, enter control acceptance criteria for automatic comparisons, and to select positive and negative interpretations based on concentration values.

- Economical open system
- Compact design
- Automatically positions each strip for reading
- User programmable test menu
- Multiple calculation modes
- Touch screen display and on-board printer
- Three strip loading capacity
- Dependable maintenance-free design
- Outputs results in 30 seconds per 12 well strip
- Multiple calculation modes
- PC connectivity via USB
- NRTL listed, CE marked, Made in the USA

 **AWARENESS[®]**
TECHNOLOGY, INC.

LumiStat[®]

Semi-automated CLIA Programmable Analyzer

Sensitivity/Detection Limit	HRP 1×10^{-18} moles, alkaline phosphatase (ALP) 1×10^{-21} moles
Linear Dynamic Range	10^6 RLU
Cross Talk	Less than 2.5×10^{-4}
Detector	Photomultiplier Tube (PMT)
Spectral Response Range	340-650 nm
Peak Wavelength	400 nm
Detection Mode	Glow luminescence
Carrier Format	Strip carrier can hold three 12-well strips or three 8-well strips
Vessel	Single, or break apart and non-break apart strips up to 12 wells long, 3 strip loading capacity. 1 – 36 well loading capacity
Display	Interactive touch screen 3.5" (8cm) LCD, color graphic display
Printer	Thermal with graphic capability, approximately 29 characters per line
Speed	Reads, calculates and prints results, in approximately 30 seconds per 12-well strip
Calculation Modes	Single Standard, Factor, multi-point calibration with point-to-point or cubic spline with curve storage, linear regressions with log, linear, and log-logit selections and cut off mode selections
Test Storage	Stores all parameters including, calculations, unit codes, ranges, interpretations, calibrator values, test names, and previous curve
Dimensions & weight	Approximately 9"x13.5"x5" (24x34x13cm) weight: 5 lbs (2.3kg)

Distribuidor Autorizado para Bolivia:



EQUIPAMIENTO MÉDICO Y DE LABORATORIO

La Paz: Av. Ormachea #125 esq. Calle 6 – Obrajes;
Telf.: (591-2) 2788493 – 72200805; cje@diamedicimport.com
Cochabamba: Av. Papa Paulo #1158 entre calles Venezuela y German Urquidi; Telf.: (591-4) 4539638- 72221300;
pje@diamedicimport.com
Santa Cruz: C. Tarija #699 esq. Saavedra – Casco Viejo;
Telf.: (591) 72200814; info@diamedicimport.com

ISO 13485:2016



314900B



Awareness Technology, Inc. USA & International
tel: +1 772 283 6540 fax: +1 772 283 8020

Awareness Technology, Inc. Europe
tel: +43 (0) 2236 892465 fax: +43 (0) 2236 892464

Awareness Technology, Inc. Asia & Africa
tel: +971 (06) 557 80 58 fax: +971 (06) 557 80 59